



SURGICAL INSTRUMENTS AND TRIALS INSTRUCTIONS FOR USE

R_x ONLY
Printed in U.S.A.

PLEASE READ THIS IN CONJUNCTION WITH SURGICAL PROCEDURE























Content	Section	Pages
Labeling Symbols Glossary		3-4
Warnings and Precautions	1	5-6
Trials and Instrument Cleaning	2	6
Trials and instrument sterilization	3	8
Re-sterilization	4	9
Declaration of contamination status	Appendix 1	11
Release Note	Appendix 2	12













Onkos Surgical Instrumentation consists of manual surgical instruments intended for use in Orthopedic surgeries.

Utilization: Before clinical use, the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the instrumentation. Surgical Protocols provide additional procedural information.

Labeling Symbols Glossary

The following is a list of symbols that may be used on Onkos Surgical medical device labeling. Refer to individual product labels for applicable symbology for each product.

Symbol		Standard
Caution (See instructions for use for operator awareness)		ISO 15223-1
Consult instructions for use or consult electronic instruction use		ISO 15223-1
Do not re-use		ISO 15223-1
Sterilized using irradiation		ISO 15223-1
Sterilized using hydrogen peroxide		N/A
Sterilized using ethylene oxide		ISO 15223-1
Sterilized using Aseptic processing techniques (Aseptic fill)		ISO 15223-1
Non-sterile		ISO 15223-1
Do not re-sterilize		ISO 15223-1
Single sterile barrier system		ISO 15223-1
Double sterile barrier system		ISO 15223-1
Use-by date		ISO 15223-1
Date of manufacture		ISO 15223-1
Legal manufacturer		ISO 15223-1
UK Responsible Person		N/A
Catalogue number		ISO 15223-1
Batch code		ISO 15223-1
Serial number		ISO 15223-1
MR Safe		ASTM F2503
MR Conditional		ASTM F2503
MR Unsafe		ASTM F2503
Do not use if package is damaged		ISO 15223-1

Symbol	Standard
Medical device 	ISO 15223-1
Quantity 	N/A
"Caution: Federal Law (USA) restricts this device to sale by on the order of a physician" 	N/A
Contains hazardous substances 	ISO 15223-1
Temperature limit 	ISO 15223-1
Keep dry 	ISO 15223-1
Keep away from sunlight 	ISO 15223-1
Patient identification 	ISO 15223-1
Date (of Implantation) 	ISO 15223-1
Health care center or doctor 	ISO 15223-1
Patient information website 	ISO 15223-1
Unique Device Identifier 	ISO 15223-1

PLEASE READ THIS IN CONJUNCTION WITH SURGICAL PROCEDURE

Section 1

WARNINGS and PRECAUTIONS



DO NOT implant the trial implants.



Onkos Surgical trial implants must not be used with products from other manufacturers.



The surgical procedure must be read prior to carrying out any surgical procedure.



For instructions in the use of the trial instruments refer to the surgical planning guide provided.



Improper use or mishandling of the components can result in damage to one or more of the components, or improper selection of implants.



The surgical instruments and trials are supplied **NON-STERILE**. Clean and sterilize before use, in accordance with the instructions provided (see Section 2 and 3).



Surgical instruments must never be placed in physiological saline solution as contact leads to pitting and corrosion.



Improper cleaning tools (e.g. overlarge brushes, metal brushes) can cause damage. **DO NOT** use scouring agents or steel wool.



Ensure detergent solutions operate within a pH range of 6.0-8.0.



Any deviation from recommended sterilization methods must be validated by the user.



The recommended sterilization method is only valid with sterilization equipment that is maintained, calibrated and operated in accordance with local regulatory guidelines.



Ensure steam quality meets acceptable standards to prevent damage and discoloration.



Do not sterilize using either Ethylene Oxide (EtO) or cold sterilization techniques.

PLEASE READ THIS IN CONJUNCTION WITH SURGICAL PROCEDURE

ADDITIONAL INFORMATION

Instruments can become damaged or may even fracture after repeated use or if used with excessive force. Instruments should be inspected prior to surgery and not used if damaged or worn excessively.

Unless otherwise specified on the instrument, markings on instruments with a measuring function have a linear accuracy of ± 0.5 mm.

Cleaning and Sterilization

All loaner Instruments and trials must be cleaned and decontaminated prior to returning to Onkos Surgical in accordance with the procedures defined in Section 2 and the form in Appendix 1 must be completed and returned with any items being sent back to Onkos Surgical.

Section 2

TRIALS and INSTRUMENTS CLEANING INSTRUCTIONS

All loaner instruments are supplied washed but are NON-STERILE and must be cleaned and sterilized before use in accordance with the instructions in Sections 2 and 3 and the form in Appendix 1 must be completed.

a. Information on methods

Listed here are validated methods for disinfection, cleaning and sterilization for user reprocessing of Onkos Surgical, Orthopedic Instruments and Trials sets. As far as possible, it is recommend that decontamination is performed using an automatic washing-disinfector utilizing thermal disinfection, followed by steam sterilization.

b. Preparation for cleaning

At the point of use, remove all excess visible soil and keep the reusable medical devices moist until cleaning. Any device capable of disassembly must be disassembled before cleaning. Ensure that thorough cleaning and rinsing is performed as soon as possible. Manual cleaning may be undertaken either by immersion or non-immersion techniques dependent upon the device construction but these instructions relate solely to a validated automated cleaning cycle.

PLEASE READ THIS IN CONJUNCTION WITH SURGICAL PROCEDURE

c. Cleaning Equipment and Agents

Water	Water quality is an important consideration in all cleaning steps. Demineralized or deionized water is recommended as this can help prevent discoloration and staining.
Detergents and Cleaning Agents	Only suitable detergents and cleaning agents within a range of pH 6.0-8.0 may be used.
Ultrasonic Cleaner	Ultrasonic cleaning is suitable for an especially thorough, mild cleaning of heavily soiled, and other difficult to access locations such as joints, crevices and channels.
Cleaning Equipment	Non-abrasive low linting clothe and general purpose cleaning brushes.

d. Product Cleaning Guidelines – Automatic

The following validated cleaning instruction is based on general product features which present challenges to the cleaning process and not on specific products. The instructions presented are based on an automated method using a washer-disinfector and have been validated. In the case of non-metals, unless otherwise instructed, reprocess following these guidelines.

Equipment and reagents:

Washer/disinfector (Medisafe SI Auto Workstation), Enzyme based detergent (3E-Zyme).

Procedure:

- i. Load the instruments into washer/disinfector ensuring that the basket is not overloaded and that hinges are open and cannulation holes are able to drain.
- ii. Run cycle (key stages):
 - Cold pre-wash: 3 minutes, < 35°C (<95°F)
 - Main wash: 13 minutes, 43°C (109°F)
 - Pre-Rinse: 2 minutes, 20°C (68°F)
 - Disinfection fill
 - Thermal disinfection: 2 minutes, 80-85°C (176-185°F)
 - Air purge dry: 2 minutes
- iii. When the cycle is complete remove the contents from the washer for inspection.
- iv. Inspect each instrument for cleanliness. If soiling is still evident repeat the cleaning cycle.

PLEASE READ THIS IN CONJUNCTION WITH SURGICAL PROCEDURE

Section 3

TRIALS and INSTRUMENT STERILISATION

IMPORTANT



All trials and instruments **MUST** be cleaned prior to sterilization (see section 2).



Partially unthread the polymeric cap of the General Impactor (approximately 2mm) while still ensuring that it fits in its specified location within the tray, to improve steam penetration during sterilization.

The following sterilization process should be used for trials and instrumentation.

The trials and instrumentation are recommended to be sterilized using pre-vacuum or porous load, high temperature steam sterilization (air removal via pulsed pre-vacuum method).

These devices must be placed in suitably wrapped porous membrane packaging for the sterilization process (i.e. central supply wrap, autoclave bags, paper/plastic pouches etc.).

The suggested cycles, based on recommended industry standards and regulatory guidance consist of the following:

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Pre-vacuum	134-137°C (273 - 278°F)	3 minutes (minimum)	30 minutes (minimum)

Method	Cycle Type	Temperature	Exposure Time	Drying Time
Steam	Porous load	132°C (270°F)	4 minutes (minimum)	30 minutes (minimum)

Onkos Surgical has validated the above recommended sterilization cycles. Other sterilization methods and cycles may also be suitable. Individuals or hospitals are advised, however, to validate whichever method is considered appropriate for their organization. Ethylene Oxide (EtO) and cold sterilization techniques are not recommended.

Ensure the polymeric cap is re-tightened on the General Impactor prior to use. Ensure all components are dry prior to use.

PLEASE READ THIS IN CONJUNCTION WITH SURGICAL PROCEDURE

Section 4


RE-STERILISATION

Orthopedic instruments and trials should be cleaned and re-sterilized in accordance with the instructions in section 2 and 3.



Manufacturer:

Onkos Surgical
77 East Halsey Road
Parsippany, NJ USA 07054
+1-973-264-5400
www.onkossurgical.com



UK Responsible Person:

IdealMed Ltd, Unit A2,
Beech House Oaklands
Office Park, Hooton
Cheshire, CH66 7NZ,
United Kingdom

APPENDIX 1

DECLARATION OF CONTAMINATION STATUS

This form is only to be used for the return of Trials and Instruments.

Local forms can be supplied as an alternative as long as all the required information is supplied.

Product Description		Product Identification		
From		Address		
Contact Name		Emergency contact number		
Contact e-mail		Have any of the items been contaminated		
Yes *	No	Don't Know	Please circle	
* State type of contamination: Blood, Body fluids, or any other hazard				
Have the items been decontaminated		Yes †	No ‡	Don't Know
† Was the process in accordance with the information supplied within section 2 of document JTS IFU06.02.23v0?		Yes	No	Please circle
IF NO please provide details of what cleaning and decontaminating process was used				
‡ Please explain why the items have not been decontaminated				
Signature of person completing the form		Date		
Print Name	Job Title			
By signing this form you are confirming that all of the information is correct and accurate to the best of your knowledge at the time of approval.				
<p>ANY CONTAMINATED INSTRUMENTS OR TRIALS <u>MUST NOT</u> BE RETURNED WITHOUT THE PRIOR AGREEMENT AND KNOWLEDGE OF ONKOS SURGICAL .</p> <p>THE INSTRUMENTS OR TRIALS MUST BE RETURNED SUITABLY PACKAGED AND IDENTIFIED.</p>				

APPENDIX 2
RELEASE NOTE

CRC Number: _____

SURGICAL INSTRUMENTS AND TRIALS: (Decontamination Certificate)

Special Instruments/ Trials included with the above system	Yes	<input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	No
--	-----	--	----

These instruments/trials were previously used in a surgical invasive procedure and were exposed to blood, body fluids or pathological samples. These items were subsequently cleaned and decontaminated prior to return to Onkos Surgical By (Insert organization name or leave blank)
In accordance with the validated process detailed in JTS IFU06.02.23v0.

After inspection all of the instruments and trials are washed.

I declare that I have taken all reasonable steps to ensure the accuracy of the above information

Signature	Date
-----------	------

Print Name

If further information is required on Onkos Surgical devices or instrumentation, please contact the Onkos Patient Solutions team:

Tel: +1-973-264-5400
Email: patientsolutions@onkossurgical.com